


CT SPHERIS 209M
Technical Specifications



We make it visible.

CT SPHERIS 209M

CT SPHERIS® 209M



Optic Design	Monofocal, spheric
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size	1.8 mm
Company Labeled A-Constant ¹	117.9
Diopter Range	From 0.0 to +32.0 D 0.0 to +10.0 D, 1.0 D increments +10.0 to +30.0 D, 0.5 D increments +30.0 to +32.0 D, 1.0 D increments
ACD	4.91
Implantation in	Bag
Injector / Cartridge Set ²	AT.Shooter A2-2000 / ACM2 <i>or</i> AT.Shooter A1-2000 / VISCOJECT™-BIO 1.8 Cartridge Set <i>or</i> VISCOJECT-BIO 1.8 Injector Set

¹Please refer to our web pages for optimized A-Constants. ²Please refer to our web pages for the most up-to-date references.



0297 CT SPHERIS 209M
ACM2 Cartridge Set



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.zeiss.com/iol
www.zeiss.com/med/contacts

EN_32_022_009811

CT SPHERIS is a registered trademark of Carl Zeiss Meditec. VISCOJECT is a trademark of Medica (0482). A6 injector is a trademark of MDJ (CE 0120). The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper.
© 2015 by Carl Zeiss Meditec AG. All copyrights reserved.

