



VISALIS® 500
Accessory catalog



We make it visible.

Table of contents	Page
I/A Tubing	3
Cataract Accessories	3
Phacoemulsification	3
Handpiece	3
Phaco Sets	4
Tips	6
Additional Accessories	6
Sleeves	7
Irrigation/Aspiration	7
Bimanual Technique	7
Coaxial Technique	8
Diathermy	10
Anterior Vitrectomy	10
Vitreoretinal Accessories	
Trocar and Infusion Systems	11
Posterior Vitrectomy	11
Endoillumination	12
Injection and Removal Systems	13
Endodiathermy	14
Miscellaneous Accessories	14
System Options	15

I/A Tubing

single-use		<p>I/A TUBING CASSETTE QUICKSET 303060-0204-000 Cassette I/A tubing "QuickSet" with drainage bag (40 pcs)</p>
------------	---	--

Cataract Accessories

Phacoemulsification







Handpiece

reusable		<p>ULITE PHACO HANDPIECE 303060-0001-000 Phaco Handpiece 4 piezoelectric, 40 KHz Tip stroke from 0 to 100 μm</p>
----------	---	--

Cataract Accessories

Phacoemulsification

Phaco Sets







single-use		PHACO SET 20G, U/S TIP 15° (10 pcs)	20G Phaco Tip Straight, 15° 20G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1000-000
single-use		PHACO SET 20G, U/S TIP 30° (10 pcs)	20G Phaco Tip Straight, 30° 20G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1001-000
single-use		PHACO SET 20G, U/S TIP 30° FLARED (10 pcs)	20G Phaco Tip Flared 30° 20G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1010-000
single-use		PHACO SET 21G, U/S TIP 15° (10 pcs)	21G Phaco Tip Biconical 15° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1002-000
single-use		PHACO SET 21G, U/S TIP 30° (10 pcs)	21G Phaco Tip Biconical 30° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1003-000
single-use		PHACO SET 21G, U/S TIP 30° FLARED (10 pcs)	21G Phaco Tip Flared 30° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1011-000

* Recommended incision size in mm

Cataract Accessories

Phacoemulsification

Phaco Sets

single-use		PHACO SET 22G, U/S TIP 15° (10 pcs) 22G Phaco Tip Biconical 15° 22G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1004-000
single-use		PHACO SET 22G, U/S TIP 30° (10 pcs) 22G Phaco Tip 30° Biconical 22G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1005-000
single-use		PHACO SET 22G, U/S TIP 30° FLARED (10 pcs) 22G Phaco Tip Flared 30° 22G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1012-000
single-use		20G SLEEVES, TEST CHAMBER SET (10 pcs) 20G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 2.6-2.8	303060-1006-000
single-use		21G SLEEVES, TEST CHAMBER SET (10 pcs) 21G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 2.2-2.5	303060-1007-000
single-use		22G SLEEVES, TEST CHAMBER SET (10 pcs) 22G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 1.7-2.1	303060-1008-000

* Recommended incision size in mm



Cataract Accessories

Phacoemulsification

Tips

reusable		19G PHACO TIP,15°	Straight 15°	303060-0050-000	RIS* 2.9-3.2	Sleeve 303060-0100-000
		19G PHACO TIP,30°	Straight 30°	303060-0051-000		
		19G PHACO TIP,45°	Straight 45°	303060-0052-000		
reusable		20G PHACO TIP,15°	Straight 15°	303060-0054-000	RIS* 2.6-2.8	Sleeve 303060-0101-000
		20G PHACO TIP,30°	Straight 30°	303060-0055-000		
		20G PHACO TIP,45°	Straight 45°	303060-0056-000		
reusable		20G PHACO TIP,BICONICAL	Biconical 15°	303060-0053-000	RIS* 2.6-2.8	Sleeve 303060-0101-000
reusable		20G PHACO TIP,20°,ANGLED	Angled 20°	303060-0062-000	RIS* 2.6-2.8	Sleeve 303060-0101-000
reusable		21G PHACO TIP,BICONICAL,15°	Biconical 15°	303060-0059-000	RIS* 2.2-2.5	Sleeve 303060-0102-000
		21G PHACO TIP,BICONICAL,30°	Biconical 30°	303060-0060-000		
reusable		22G PHACO TIP,BICONICAL,15°	Biconical 15°	303060-0063-000	RIS* 1.7-2.1	Sleeve 303060-0103-000
		22G PHACO TIP,BICONICAL,30°	Biconical 30°	303060-0065-000		
reusable		22G PHACO TIP,30°, SLEEVELESS BIMANUAL	Straight 30°	303060-0057-000	RIS* 1.0-1.2	
reusable		22G PROGR. PHACO TIP, SLEEVELESS BIMANUAL	Flared 30°	303060-0058-000	RIS* 1.0-1.2	
reusable		PARS PLANA PHACO TIP	20G (Posterior Fragmentation)	303060-0064-000		

Additional Accessories

reusable		TEST CHAMBER (2 pcs)		303060-0202-000
reusable		PHACO TIP WRENCH		303060-0121-000

* Recommended incision size in mm

Cataract Accessories

Phacoemulsification

Sleeves

reusable		19G SILICONE SLEEVE (GREY) 2.9-3.2MM (3 pcs)	RIS* 2.9-3.2	303060-0100-000
reusable		20G SILICONE SLEEVE (BLUE) 2.6-2.8MM (3 pcs)	RIS* 2.6-2.8	303060-0101-000
reusable		21G SILICONE SLEEVE (YELLOW) 2.2-2.5MM (3 pcs)	RIS* 2.2-2.5	303060-0102-000
reusable		22G SILICONE SLEEVE (GREEN) 1.7-2.1MM (3 pcs)	RIS* 1.7-2.1	303060-0103-000

Irrigation/Aspiration

Bimanual Technique

single-use		IRR & ASP HANDPIECES BIMANUAL (10 pcs) Irrigation Side Ports & Aspiration Bimanual Handpieces Aspiration: 21G, 0.35 mm Irrigation: 21G, 2 Side Ports 0.45 mm		303060-0165-000
single-use		IRR & ASP HANDPIECES BIMANUAL FRONT (10 pcs) Irrigation Front Port & Aspiration Bimanual Handpieces Aspiration: 21G, 0.35 mm Irrigation: 21G, 1 Front Port 0.5 mm		303060-0166-000
reusable		IRRIGATION HANDPIECE, BIMANUAL 21G, 0.6 mm		303060-0155-000
reusable		ASPIRATION HANDPIECE, BIMANUAL 21G, 0.3 mm		303060-0156-000

* Recommended incision size in mm

Cataract Accessories

Irrigation/Aspiration

Coaxial Technique

single-use		I/A HANDPIECE 20G STRAIGHT 20G I/A Handpiece (10 pcs), Coaxial, Straight	RIS* 2.6-2.8	303060-0141-000
single-use		I/A HANDPIECE 20G ANGLED 45° 20G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 2.6-2.8	303060-0142-000
single-use		I/A HANDPIECE 21G ANGLED 45° 21G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 2.2-2.5	303060-0143-000
single-use		I/A HANDPIECE 21G CURVED 21G I/A Handpiece (10 pcs), Coaxial, Curved 40°	RIS* 2.2-2.5	303060-0146-000
single-use		I/A HANDPIECE 22G ANGLED 45° 22G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 1.7-2.1	303060-0144-000
single-use		I/A HANDPIECE 22G CURVED 22G I/A Handpiece (10 pcs), Coaxial, Curved 40°	RIS* 1.7-2.1	303060-0145-000
reusable		I/A HANDPIECE, COAXIAL, REUSABLE I/A Handpiece, Coaxial		303060-0140-000

* Recommended incision size in mm

Cataract Accessories

Irrigation/Aspiration




Coaxial Technique

reusable		1/A CANNULA, METAL SLEEVE, STRAIGHT, Ø1.6 MM		303060-0150-000
		1/A CANNULA, METAL SLEEVE, ANGLED, Ø1.6 MM		303060-0151-000
		1/A CANNULA, METAL SLEEVE, CURVED, Ø1.6 MM		303060-0152-000
		1/A CANNULA, MET. SLEEVE, 12 O'CLOCK, Ø1.6 MM		303060-0153-000
		1/A CANNULA, METAL SLEEVE, 90°, Ø1.6 MM		303060-0154-000
reusable		19G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.9-3.2	303060-0157-000
		19G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0158-000
reusable		20G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.6-2.8	303060-0159-000
		20G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0160-000
reusable		21G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.2-2.5	303060-0161-000
		21G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0162-000
reusable		22G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 1.7-2.1	303060-0163-000
		22G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0164-000

* Recommended incision size in mm

Cataract Accessories

Diathermy

reusable		DIATHERMY PENCIL ERASER	303060-0252-000
reusable		DIATHERMY FORCEPS	303060-0251-000
reusable		DIATHERMY BIPOLAR CABLE	303060-0250-000

Anterior Vitrectomy

single-use		20G A-VIT PROBE WITH SLEEVE (5 pcs) 20G Anterior Vitrectomy Probe Silicone Sleeve	303060-0301-000
reusable		20G A-VIT PROBE W/SLEEVE, REUSABLE 20G Anterior Vitrectomy Probe Silicone Sleeve	303060-0300-000







Vitreoretinal Accessories

Trocar and Infusion Systems

single-use		20G SCLERAL INFUSION CANNULA (5 pcs)	303060-0503-000
single-use		23G TROCAR KIT ONE-STEP (10 pcs)	303060-0405-000
single-use		23G TROCAR KIT ONE-STEP OIL INFUSION (10 pcs)	303060-0406-000
single-use		25G TROCAR KIT (5 pcs) One-Step	303060-0409-000











Vitreoretinal Accessories

Posterior Vitrectomy

single-use 	 20G POSTERIOR VITRECTOMY PROBE (5 pcs) 303060-0403-000 High Speed
single-use 	 23G POSTERIOR VITRECTOMY PROBE (5 pcs) 303060-0407-000 High Speed
single-use 	 25G POSTERIOR VITRECTOMY PROBE (5 pcs) 303060-0408-000 High Speed





Vitreoretinal Accessories

Endoillumination

single-use		 20G ENDO-ILLUMINATION PROBE (5 pcs) 303060-0600-000
		 20G WIDE-ANGLE ENDO-ILLUMINATION PROBE (5 pcs) 303060-0601-000
single-use		 20G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs) 303060-0602-000
		 20G WIDE-ANGLE ENDO-ILLUMIN. MULTIPOINT (5 pcs) 303060-0603-000
single-use		 23G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs) 303060-0604-000
		 23G WIDE-ANGLE ENDO-ILLUMIN. MULTIPOINT (5 pcs) 303060-0605-000
single-use		25G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs) 303060-0606-000

Vitreoretinal Accessories

Injection and Removal Systems

single-use		ADAPTER AIR / CONTROLLED IRRIGATION (5 pcs)	303060-0506-000
single-use		AIR INJECTION TUBE WITH FILTER (5 pcs)	303060-0508-000
single-use		SILICONE-OIL INJECTION TUBE (5 pcs)	303060-0509-000
single-use		VISCOUS FLUID REMOVAL KIT (10 pcs)	303060-0507-000

Vitreoretinal Accessories


Endodiathermy

single-use		20G ENDODIATHERMY PROBE (5 pcs)	303060-0262-000
single-use		23G ENDODIATHERMY PROBE (5 pcs)	303060-0263-000
single-use		25G ENDODIATHERMY PROBE (5 pcs)	303060-0264-000
reusable		DIATHERMY BIPOLAR CABLE	303060-0250-000





Miscellaneous Accessories

single-use		SCREEN DRAPE VISALIS 500 (8 pcs)	303060-0352-000
single-use		REMOTE CONTROL COVER VISALIS (10 pcs)	303060-0353-000
single-use		TRAY COVER VISALIS 500 (2 pcs)	303060-0354-000

Miscellaneous Accessories

reusable	 A rectangular blue and green sterilization tray with gold-colored contacts on the top surface.	STERILIZATION TRAY 26,5 x 16,3 x 4,6 cm	303060-0350-000
reusable	 A white, cylindrical touchscreen pen with a small circular button at the top.	TOUCHSCREEN PEN (5 pcs)	301640-0020-000

System Options

	 A white, rectangular console tray with a metal handle on the right side and a small latch on the left.	CONSOLE TRAY with manually adjustable height	303041-9004-000
	 A white, rectangular video overlay box with a blue logo on the front panel.	VIDEO OVERLAY BOX	303041-9007-000
	 A white, rectangular IV-pole footswitch with two blue buttons and a cable attached to the right side.	IV-POLE FOOTSWITCH	303041-9006-000
	 A white, oval-shaped remote control with a cable attached to the top and a blue logo on the front.	REMOTE CONTROL	303041-9008-000

Your local contact:**Argentina**

Carl Zeiss Argentina S.A.
Calle Nahuel Huapi 4015 / 25
C1430 BCO Buenos Aires
Argentina
Phone: +54 11 45 45 66 61
bruzzi@zeiss.com.ar

Australia

Carl Zeiss Pty Ltd
Tenancy Office 4, Level 1
40-52 Talavera Road
North Ryde NSW 2113
Australia
Phone: +61 2 9020 1333
med@zeiss.com

Austria

Carl Zeiss GmbH
Laxenburger Str. 2
1100 Vienna
Austria
Phone: +43 1 79 51 80
austria@zeiss.org

Belgium

Carl Zeiss NV-SA
Ikarooslaan 49
1930 Zaventem
Belgium
Phone: +32 2 719 39 11
info@zeiss.be

Brazil

Carl Zeiss do Brasil Ltda.
Av. Nações Unidas, 21711
CEP04795-100 São Paulo
Brazil
Phone: +55 11 5693 5521
medbrasil@zeiss.org

Canada

Carl Zeiss Canada Ltd.
45 Valleybrook Drive
Toronto, ON M3B 2S6
Canada
Phone: +1 800 387 8037
micro@zeiss.com

China

Carl Zeiss Shanghai Co. Ltd.
1/f., Ke Yuan Building
11 Ri Yin Nan Road
Waigaoqiao Free Trade Zone
2005 Yang Gao Bei Road
Shanghai 200131
China
Phone: +86 21 5048 17 17
sro@zeiss.com.cn

Czech Republic

Carl Zeiss spol. s r.o.
Radlická 14/3201
150 00 Prague 5
Czech Republic
Phone: +420 233 101 221
zeiss@zeiss.cz

France

Carl Zeiss Meditec France SAS
60, route de Sartrouville
78230 Le Pecq
France
Phone: +33 1 34 80 21 00
med@zeiss.fr

Germany

Carl Zeiss Meditec VG mbH
Carl-Zeiss-Strasse 22
73447 Oberkochen
Germany
Phone: +49 7364 20 6000
vertrieb@meditec.zeiss.com
Surgical Ophthalmology:
Phone: +49 800 470 50 30
iol.order@meditec.zeiss.com

Hong Kong

Carl Zeiss Far East Co. Ltd.
Units 11-12, 25/F
Tower 2, Ever Gain Plaza
No. 88 Container Port Road
Kwai Chung
Hong Kong
Phone: +852 2332 0402
cffe@zeiss.com.hk

India

Carl Zeiss India Pvt. Ltd.
Plot No.3, Jigani Link Road
Bommasandra Industrial Area
Bangalore - 560 099
India
Phone: +91 80 4343 8000
info@zeiss.co.in

Italy

Carl Zeiss S.p.A.
Viale delle Industrie 20
20020 Arese (Milan)
Italy
Phone: +39 02 93773 1
infomed@zeiss.it

Japan

Carl Zeiss Meditec Japan Co. Ltd.
Shinjuku Ku
Tokyo 160-0003
22 Honchio-Cho
Japan
Ophthalmic instruments:
Phone: +81 3 33 55 0331
medsales@zeiss.co.jp
Surgical instruments:
Phone: +81 3 33 55 0341
cmskoho@zeiss.co.jp

Malaysia

Carl Zeiss Sdn Bhd.
Lot2, Jalan 243/51 A
46100 Petaling Jaya
Selangor Darul Ehsan
Malaysia
Phone: +60 3 7877 50 58
malaysia@zeiss.com.sg

Mexico

Carl Zeiss de México S.A. de C.V.
Avenida Miguel Angel de Quevedo
496
04010 Mexico City
Mexico
Phone: +52 55 59 99 0200
cz-mexico@zeiss.org

Netherlands

Carl Zeiss B.V.
Trapezium 300
Postbus 310
3364 DL Sliedrecht
Netherlands
Phone: +31 184 43 34 00
info@zeiss.nl

New Zealand

Carl Zeiss (N.Z.) Ltd.
15B Paramount Drive
P.O. Box 121 - 1001
Henderson, Auckland 0650
New Zealand
Phone: +64 9 838 5626
med@zeiss.com

Poland

Carl Zeiss sp. z o.o.
ul. Lopuzanska 32
02-220 Warsaw
Poland
Phone: +48 22 858 2343
medycyna@zeiss.pl

Singapore

Carl Zeiss Ptd. Ltd.
50 Kaki Bukit Place
Singapore 415926
Singapore
Phone: +65 6741 9600
info@zeiss.com.sg

South Africa

Carl Zeiss (Pty.) Ltd.
363 Oak Avenue
Ferndale
Randburg 2194
South Africa
Phone: +27 11 886 9510
info@zeiss.co.za

South Korea

Carl Zeiss Co. Ltd.
Seoul 121-828
Mapo-gu
141-1, Sangsu-dong
2F, BR Elitel Bldg.
South Korea
Phone: +82 2 3140 2600
korea@zeiss.co.kr

Spain

Carl Zeiss Meditec Iberia S.A.U.
Ronda de Poniente, 15
Tres Cantos
28760 Madrid
Spain
Phone: +34 91 203 37 00
info@zeiss.es

Sweden

Carl Zeiss AB
Tegeluddsvaegen 76
10254 Stockholm
Sweden
Phone: +46 84 59 25 00
info@zeiss.se

Switzerland

Carl Zeiss AG
Feldbachstrasse 81
8714 Feldbach
Switzerland
Phone: +41 55 254 7534
info@zeiss.ch

Thailand

Carl Zeiss Thailand
Floor 8, Thosapol Land Building 2
230 Ratchadapisek Road
Huaykwang, Bangkok 10310
Thailand
Phone: +66 2 2 74 06 43
thailand@zeiss.com.sg

United Kingdom

Carl Zeiss Ltd.
509 Coldhams lane
Cambridge
CAMBS
CB1 3JS,
United Kingdom
Phone: +44 1707 871200
info@zeiss.co.uk

United States of America

Carl Zeiss Meditec, Inc.
5160 Hacienda Drive
Dublin, CA 94568
USA
Phone: +1 925 557 4100
info@meditec.zeiss.com

EN_30_010_01131 Printed in Germany CZ-VIII/2012
The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information.
Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper.
VISALUS is a registered trademark of Carl Zeiss. © 2011 by Carl Zeiss Meditec AG. All copyrights reserved.

